

Vilistus Mind Mirror 6

Self-Training Kit



Vilistus-4 Digital Sampling Unit (DSU)

Use the DSU's four instrument ports for EEG and physiology monitoring and/or the optional Manual Trigger Button. Choose Bluetooth or USB connectivity to the computer. Takes two 1.5-volt AA lithium batteries, not supplied. Two-year fair use warranty.



Dual-Channel EEG Sensor

Snap in four reusable flat blue disk electrodes for two active and two passive sensor heads and one pre-gelled disposable electrode for easy ground wire placement. (See electrodes here and below.) Shielded wires reduce electrical artifact. Plugs into a single DSU port. One-year fair use warranty.



Reusable Electrodes

Flat metal disk electrodes made by Multibiosensors. 50 supplied. Good for 3-8 EEG uses by the same person, when gently cleaned with tissues.



Disposable Electrodes. Covidien pre-gelled electrodes for easy EEG ground wire placement on the forehead, neck, chest or shoulders. 50 supplied. For single use.



NuPrep Abrasive Gel

A small amount is used on a Q-tip to slightly abrade and clean the scalp before applying the electrodes. 4-ounce tube.



Ten20 Conductive Paste

Apply to EEG electrodes to hold them in place on the scalp and transmit electrical signals. This stiff paste contains an excellent balance of adhesiveness and conductivity. 4-ounce tube.



Hard-shell Carry Case

This durable 12x12.5x3-inch carry case, made of recyclable polypropylene, contains a foam filling that keeps equipment secure while stored or on the move. Catches close tightly. Comfortable handle. Foam can be easily cut to custom-fit your equipment.